



Star wedges

Never falling off the wedge

Star wedges has a rhombus cutout in the tip allowing a strong resistance at lingual to avoid falling off during operation.

Star wedges with flexible material to avoid trauma to the gingival of the patients.

Hour-glass shape for easy placement

Large tab with a hole for facilitating placement and removal with hand or a pin tweezers.



4 wedge size:
Extra small, small, medium, large

Rhombus cutout in the tip allows easy insertion and causes less trauma to patients

Star wedges can squeeze and expand back after insertion

Better hold to avoid falling off during operation

